


Physics and Chemistry Lab Supplies





We offer a variety of magnetic compasses, from small plotting compasses up to a 45mm compass with a cover and locking mechanism. Our compasses include sturdy aluminum cases and magnetic needles with the North pole marked.

CATALOG NO.	DESCRIPTION	SIZE
MLS-C1SG12	Magnetic Compasses with Aluminum case and glass top	12 mm
MLS-C1SG16	Magnetic Compasses with Aluminum case and glass top	16 mm
MLS-CBSG18	Magnetic Compasses with Aluminum ring, glass top and bottom	18 mm
MLS-C1SG20	Magnetic Compasses with Aluminum case and glass top	20 mm
MLS-C1SG25	Magnetic Compasses with Aluminum case and glass top	25 mm
MLS-CWR025	Magnetic Compasses with Aluminum case, glass top with ring	25 mm
MLS-CWR030	Magnetic Compasses with Aluminum case, glass top with ring	30 mm
MLS-CPL035	Magnetic Compasses with Plastic case and plastic top	35 mm
MLS-CWR045-W	Magnetic Compasses with Aluminum case, glass top with ring	45 mm
MLS-CWR045-B	Magnetic Compasses with Aluminum case, glass top with ring	45 mm
MLS-CWL045	Magnetic Compasses with Aluminum case, glass top and lock	45 mm

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-TTHW01	Wooden Test Tube Rack. Six 25mm holes, six drying pins, with a 1/2" thick base.

CATALOG NO.	DESCRIPTION	DIMENSIONS
MLS-OFLSET	Aluminum overflow can and catch bucket.	Can: 4" I.D. X 5" H, CAP. 800ML, Spout: 0.5" I.D. X 1.5" L Bucket: 4" I.D. X 3" H
The can and bucket are also sold separately.		



CATALOG NO.	DESCRIPTION
MLS-ACLP01	Alligator Clips; 2" long, plain, chrome plated, with excellent conductivity. Available in packs of 10 or 100.

CATALOG NO.	DESCRIPTION
MLS-LBJK44	Laboratory Jacks with 4" x 4" platform
MLS-LBJK66	Laboratory Jacks with 6" x 6" platform
MLS-LBJK88	Laboratory Jacks with 8" x 8" platform
MLS-LBJK12	Laboratory Jacks with 12" x 12" platform
Laboratory support jacks are constructed of painted aluminum and stainless steel. The scissor jack design enables height adjustment with a simple turn of the knob.	



CATALOG NO.	DESCRIPTION
MLS-BUBR01	Burette Brush; Nylon bristles mounted in twisted heavy gauge galvanized wire. Includes a radial top.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-CCK001	Create-A-Circuit Kit. It comes complete with a battery holder (D size battery not included), a knife switch, two miniature light bulbs, two lamp holders, and four alligator clip connecting cords. Includes activity guide.

CATALOG NO.	DESCRIPTION
MLS-CLBUR1N	Burette Clamp, Single
MLS-CLBUR2N	Burette Clamp, Double
These clamps accept burettes up to 35mm O.D., and can be used with support rods up to 1/2" O.D.	



CATALOG NO.	DESCRIPTION
MLS-CLEX01	2-Prong Extension Clamp with Steel Rod. Features a 2-prong die-cast clamp attached to a chrome-plated steel rod. Jaws are lined with cork. Overall length is 8".

CATALOG NO.	DESCRIPTION
MLS-CLLB02	Lincoln-Type Burette Clamp. It can be used with rods up to 3/4" in diameter. The distance between the two burettes is 7-1/4" (18.4cm).



CATALOG NO.	DESCRIPTION
MLS-CLMKE1	Knife Edge Lever Clamp

CATALOG NO.	DESCRIPTION
MLS-CLMKE1	Plated steel double clamp. Accepts rods up to 3/4" diameter (19 mm).



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-CLMP01	Mohr's Pinchcock Tubing Clamp, for 10mm Diameter Tubing.
MLS-CLMP02	Mohr's Pinchcock Tubing Clamp, for 13mm Diameter Tubing.
MLS-CLMP03	Mohr's Pinchcock Tubing Clamp, for 16mm Diameter Tubing.

CATALOG NO.	DESCRIPTION
MLS-CLMPC2	C-Clamp; 2" jaw opening.
MLS-CLMPC3	C-Clamp; 3" jaw opening.
MLS-CLMPC4	C-Clamp; 4" jaw opening.
MLS-CLMPC5	C-Clamp; 5" jaw opening.
MLS-CLMPC6	C-Clamp; 6" jaw opening.
MLS-CLMPC7	C-Clamp; 7" jaw opening.



CATALOG NO.	DESCRIPTION
MLS-CLPEN1	Pendulum Clamp

CATALOG NO.	DESCRIPTION
MLS-CLRM01	Double Walled Calorimeter



CATALOG NO.	DESCRIPTION
MLS-CLPL01	Pulley with Clamp. The sheave is 50mm in diameter. Clamp can be mounted on surfaces up to 35mm thick.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION	DIMENSIONS
MLS-CBM050	Ceramic Magnet bars	50mm x 10mm x 7mm
MLS-CBM070	Ceramic Magnet bars	70mm x 15mm x 8mm
MLS-CM7912	Ceramic Magnet disc	25mm x 5mm
MLS-CM7702	Ceramic Rectangular Magnet	27mm x 8mm x 3mm
MLS-CM7907	Ceramic Magnet ring	38mm x 8mm
MLS-CM7907-RB	Ceramic Magnet rings; pack of two. One red and one blue	38mm x 8mm



CATALOG NO.	DESCRIPTION
MLS-CBST03	Cork Borers with plastic handles.
MLS-3452	Cork Borers with solid Brass handles, set of 6
MLS-CBST09	Cork Borers with solid Brass handles, set of 9
MLS-CBST12	Cork Borers with solid Brass handles, set of 12

CATALOG NO.	DESCRIPTION
MLS-LMPR02	Miniature Lamp Holder, Fahnestock Clips



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-KSW001-K	Knife switch with plastic base. suitable for low voltage experiments (less than 3.6V). Includes screw connecting posts

CATALOG NO.	DESCRIPTION
MLS-LMPP01	Miniature Lamp Holder, Screw Type



CATALOG NO.	DESCRIPTION
MLS-LMPB01	Lamp Board with 1 Lamp
MLS-LMPB03	Lamp Board with 3 Lamp
MLS-LMPB04	Lamp Board with 4 Lamp
MLS-LMPB05	Lamp Board with 5 Lamp

CATALOG NO.	DESCRIPTION
MLS-LMP015	Miniature Lamp Bulbs, 1.5V, 0.3A, pack of 10
MLS-LMP025	Miniature Lamp Bulbs, 2.5V, 0.3A, pack of 10
MLS-LMP038	Miniature Lamp Bulbs, 3.8V, 0.3A, pack of 10
MLS-LMP062	Miniature Lamp Bulbs, 6.2V, 0.5A, pack of 10



CATALOG NO.	DESCRIPTION
MLS-KSE001	Economy Switch. Simple switch uses fahnestock clips for easy, quick connections.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-BTHD01	D Cell Battery Holder and Joiner



CATALOG NO.	DESCRIPTION
MLS-BTHP01-F	D Cell Battery Holder with Fahnestock Clips



CATALOG NO.	DESCRIPTION
MLS-KSW001	Knife Switches, Screw Posts, Single pole, Single throw
MLS-KSW002	Knife Switches, Screw Posts, Single pole, Double throw

CATALOG NO.	DESCRIPTION
MLS-INPL01	Inclined Plane with Pulley
MLS-INPL01-R	Roller only for Inclined Pulley
This versatile and tough steel inclined plane includes a graduated protractor and an adjustable pulley. Incline is 4" wide x 23" long. Includes a 14 oz.	



CATALOG NO.	DESCRIPTION
MLS-1950	Motor Generator Unit

CATALOG NO.	DESCRIPTION
MLS-1117	Light Box with accessories (Optical kit Accessories)



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-656	Linear Expansion Apparatus.

CATALOG NO.	DESCRIPTION
MLS-1114-5	Ray Box



CATALOG NO.	DESCRIPTION
MLS-1478	Whimshurst Generator.

CATALOG NO.	DESCRIPTION
MLS-1783-3	Resistance Box



CATALOG NO.	DESCRIPTION
MLS-605	Incline Plane complete with arc roller metal block
The wooden inclined board is 60cm long and 10cm wide. It can be set at any angle between 0 degrees and 45 degrees.	

CATALOG NO.	DESCRIPTION
MLS-1129	Newton Color Disc with Inst. Sheet (L-12) Black Lead Free Powder Coated Base



PHYSICS AND CHEMISTRY LAB SUPPLIES




CATALOG NO.	DESCRIPTION
MLS-1123	Newtons Ring Appts. with Inst. Sheet

CATALOG NO.	DESCRIPTION
MLS-1433	Oersteds Apparatus



CATALOG NO.	DESCRIPTION
MLS-480	Wall Mount Aneroid Barometer

CATALOG NO.	DESCRIPTION
MLS-3361-2	Burner Stand 4 ring x 9 leg



CATALOG NO.	DESCRIPTION
MLS-659	Ball Ring Double Chain Type 25mm Dia Bal

CATALOG NO.	DESCRIPTION
MLS-3201-10	Centrifuge Machine Hand 4 Tube



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-1388	Electromagnet Dissectible 100mm

CATALOG NO.	DESCRIPTION
MLS-607-1	Friction Block and Surface Set



CATALOG NO.	DESCRIPTION
MLS-80-200	Slotted Weight Brass 200gm
MLS-80-50	Slotted Weight Brass 50gm

CATALOG NO.	DESCRIPTION
MLS-3372-3	Support Ring with Clamp 4 (L-12)
MLS-3372-2	Support Ring with Clamp 3 inch O.D



CATALOG NO.	DESCRIPTION
MLS-37	Spring Balance Transparent 1000gm

CATALOG NO.	DESCRIPTION
MLS-1298-2	Tuning Fork Set of 8 in Thermocole Box



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-257	Balls Metal 25mm Set of 6 Pendulum Bob

CATALOG NO.	DESCRIPTION
MLS-1092	Magnifier Hand Lens Metal Mount 62mm Magnifier



CATALOG NO.	DESCRIPTION
MLS-3826-R	Pipette Filler rubber bulb 60ml 3-valve Red

CATALOG NO.	DESCRIPTION
MLS-1011	Prism Glass Equilateral 50mm



CATALOG NO.	DESCRIPTION
MLS-3573-3	Spatula with Wooden Handle SS 4 inch S x 19 mm

CATALOG NO.	DESCRIPTION
MLS-220	Tuning Fork Set of 8 in Thermocole Box



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-210	Micrometer; Stainless Steel screw (rod) accurately machined. Brass body, dull nickel plated finish. 25 mm x 1 mm pitch.



CATALOG NO.	DESCRIPTION
MLS-6450	Pipette Support. Horizontal for 12 pipettes. Polypropylene.



CATALOG NO.	DESCRIPTION
MLS-3826-R	Pipette Support. Round, for 28 pipettes.

CATALOG NO.	DESCRIPTION
MLS-B-161	Slide Box.



CATALOG NO.	DESCRIPTION
MLS-3014-2	Alcohol Burner. Flint glass, 120 ml capacity.

CATALOG NO.	DESCRIPTION
MLS-191/1	Calipers Vernier Steel NP 125mm in pouch



PHYSICS AND CHEMISTRY LAB SUPPLIES





CATALOG NO.	DESCRIPTION	BASE	ROD
MLS-3353-611	Support Stand With Rod	6" x 11"	1/2" x 36"
MLS-3353-64	Support Stand With Rod	4" x 6"	5/16" x 18"
MLS-3353-69	Support Stand With Rod	6" x 9"	1/2" x 24"
Support Stand made from cast iron, with an enamel finish that is acid resistant. Includes a zinc plated steel rod that screws into the base support. Available in different sizes			

CATALOG NO.	DESCRIPTION
MLS-3373-1	Wire Support Ring with Clamp; 2dia, 50mm
MLS-3373-2	Wire Support Ring with Clamp; 3dia, 75mm
MLS-3373-3	Wire Support Ring with Clamp; 4dia, 100mm
MLS-3373-4	Wire Support Ring with Clamp; 5dia, 125mm



CATALOG NO.	DESCRIPTION
MLS-247	Aluminum ball with Hook 25mm
MLS-250	Steel ball with Hook 25mm

CATALOG NO.	DESCRIPTION
MLS-3014-2	Beaker Tong. Made from 1/4" thick steel rod.



CATALOG NO.	DESCRIPTION
MLS-3513	Crucible Tong Iron N.P. 200 MM

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-3511	Utility Tongs 210 mm

CATALOG NO.	DESCRIPTION
MLS-3580	Lab Scoops. Stainless steel, 6" long.



CATALOG NO.	DESCRIPTION
MLS-3545	Test Tube Holder Big with Finger Grip

CATALOG NO.	DESCRIPTION	SIZE
MLS-3568-150	Spatula Stainless Steel	150mm x 18mm
MLS-3568-200	Steel ball with Hook 25mm	200mm x 18mm



CATALOG NO.	DESCRIPTION
MLS-3470	Deflagration Spoon

CATALOG NO.	DESCRIPTION	DIMENSIONS
MLS-13742	Magnet - Bar	75mm x 11mm x 6mm



CATALOG NO.	DESCRIPTION
MLS-1394	Magnetic Compass; Small, 12.5 mm diameter.

PHYSICS AND CHEMISTRY LAB SUPPLIES




CATALOG NO.	DESCRIPTION
MLS-HTCP01	Hand Tally Counter

CATALOG NO.	DESCRIPTION
MLS-HSCOPI	Screw Compressor Tubing Clamps; 3/4"
MLS-HSCOP2	Screw Compressor Tubing Clamps; 1"



CATALOG NO.	DESCRIPTION
MLS-HKCL01	Hook Collar Clamps; for 1/2" dia. rods
MLS-HKCL02	Hook Collar Clamps; for 3/4" dia. rods

CATALOG NO.	DESCRIPTION
MLS-GSLT01	Gas Lighter
MLS-GSFLT1-PK/5	Replacement flints; pack of 5



CATALOG NO.	DESCRIPTION
MLS-FWB500	Fractional Weight Set with rider and forceps
MLS-FWS500	Fractional Weight Set with rider
MLS-AWR005	Riders, 5G; pack of 3

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-FRSS14	Forceps, Stainless Steel, Blunt, 5"
MLS-FOBL05	Forceps, Stainless Steel, Economy Blunt, 5"
MLS-FOBLC5	Forceps, Stainless Steel, Angular Blunt, 5"
MLS-FOMF05	Forceps, Stainless Steel, Sharp, 5"
MLS-FOBL45	Forceps, Stainless Steel, Economy Blunt, 4.5"
MLS-BNCT05	Forceps, Bone cutting, 5"
MLS-BNCT07	Forceps, Bone cutting, 7"
MLS-FOPN04	Forceps, Dissecting, 4"

CATALOG NO.	DESCRIPTION
MLS-FRSG12	Friction Rod, solid glass, 12" long
MLS-FRHG12	Friction Rod, hollow glass, 12" long
MLS-FRSP12	Friction Rod, plastic, 12" long
MLS-FRHR12	Friction Rod, Hard Rubber, 12" long
MLS-FRHR06	Friction Rod, Hard Rubber, 6" long
MLS-RPSL30	Friction Pad, silk, 12" x 12"
MLS-RPCT30	Friction Pad, Cotton, 12" x 12"
MLS-RPWL30	Friction Pad, Wool, 12" x 12"



CATALOG NO.	DESCRIPTION
MLS-FSCP01-PK/10	Fahnestock Clips. Supplied as a package of 10 clips.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-WDSET1	Wave Demonstrator Set. Set includes two helical springs. The long, narrow spring measures 3/4" x 72". The larger spring measures 3" x 4"
The springs are sold separately as well.	

CATALOG NO.	DESCRIPTION
MLS-FAP075	Equilateral Refraction Prism, Acrylic, measuring 75mm long by 9mm thick.



CATALOG NO.	DESCRIPTION
MLS-EMGT01	Large U-Shaped Electromagnet. Operating on 1.5 to 6 V DC, this electromagnet lifts up to 3.5kg.

CATALOG NO.	DESCRIPTION
MLS-EGA001	Boyle's Law Apparatus. Includes activity guide.



CATALOG NO.	DESCRIPTION
MLS-ECLJ03	3-Prong Extension Clamp with Steel Rod. Prongs are covered with rubber sleeves. Overall length is 10".

CATALOG NO.	DESCRIPTION
MLS-DCM015	Miniature DC Motor




PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-DPS001	Dissecting Tray, 11.5" X 7.5" X 2"
MLS-DPS002	Dissecting Tray, 12" X 10" X 2"
MLS-DPS003	Dissecting Tray, 12.75" X 9.25" X 2"
MLS-DPS004	Dissecting Tray, 16" X 12.75" X 2"

CATALOG NO.	DESCRIPTION
MLS-DNAM01-K4	DNA Model Kit with four models and CD
MLS-DNAM02-CD	CD only for DNA model
MLS-DNAM03-U	DNA model. One Unassembled model



CATALOG NO.	DESCRIPTION
MLS-DNAKIT	DNA Manipulatives Kit. Includes 78 pieces of specially designed plastic manipulatives. Also included are a manual and instructional videotape

CATALOG NO.	DESCRIPTION
MLS-DFG003	Demo Diffraction Grating. 20mm x 10mm gratings with 100, 300, and 600 lines per mm



CATALOG NO.	DESCRIPTION
MLS-DCSET10	Density Cube Set of 10. This set of 1" cubes includes one each of acrylic, oak, pine, rosewood, steel, aluminum, copper, brass, nylon, and PVC cubes packaged in a wooden storage box.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-CTSP09	Crucible Tongs; Plated Steel
MLS-CTOS09	Crucible Tongs; Oxidized Steel
MLS-CTBR09	Crucible Tongs; Brass
MLS-CTSS09	Crucible Tongs; Stainless steel

CATALOG NO.	DESCRIPTION
MLS-COBR3-U	Burette Clamp with Boss Head, Uncoated Jaws



CATALOG NO.	DESCRIPTION
MLS-LSCSET	Laboratory Scoop and handle set.
MLS-LSCSS6	Laboratory Scoop, Stainless Steel. Length 6-1/2"
MLS-LSHND1	5" long hardwood Scoop handle

CATALOG NO.	DESCRIPTION
MLS-LSTA38	38 mm Demonstration Lens Set-Acrylic
MLS-LSTA50	50 mm Demonstration Lens Set-Acrylic
MLS-LSTA75	75 mm Demonstration Lens Set-Acrylic



CATALOG NO.	DESCRIPTION
MLS-LYJ001	Leyden Jar. Two aluminum jars (one with an attached aluminum electrode) and one plastic jar are used to collect the charge. Assembled device is 6" tall and 3-1/4" in diameter. Includes activity guide.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-MAHB03	Brain Model, 3-Part. Includes a sturdy base and a detailed key.

CATALOG NO.	DESCRIPTION
MLS-MAM001	DC Ammeters; 0-50MA/0-500MA / 0-5A
MLS-MAM002	DC Ammeters; 0-1A/0-5A
MLS-MAM003	DC Ammeters; 0-1A
MLS-MAM004	DC Ammeters; 0-3A
MLS-MAM005	DC Ammeters; 0-5A
MLS-MAM006	DC Ammeters; 0-10A



CATALOG NO.	DESCRIPTION
MLS-CPT015	Clay Pipe Triangles; 1 1/2" side length
MLS-CPT020	Clay Pipe Triangles; 2" side length
MLS-CPT025	Clay Pipe Triangles; 2 1/2" side length
MLS-CPT030	Clay Pipe Triangles; 3" side length

CATALOG NO.	DESCRIPTION
MLS-CMSET6	Floating Magnets Set. Set includes six brightly colored magnetic rings, plastic stand, and cardboard storage box.



CATALOG NO.	DESCRIPTION
MLS-CMSET1	Ceramic Disc Magnets. Magnets are 25mm in diameter.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-MBMS01	Monocotyledon Stem Model

CATALOG NO.	DESCRIPTION
MLS-PUDS01	Demonstration Single Pulley
MLS-PUDD01	Demonstration Double Pulley
MLS-PUDT01	Demonstration Triple Pulley



CATALOG NO.	DESCRIPTION
MLS-PTAN01	Single Tandem Pulley; 50mm Sheave
MLS-PTAN02	Double Tandem Pulley; 40 and 50 mm Sheaves
MLS-PTAN03	Triple Tandem Pulley; 30 mm, 40 mm and 50 mm Sheaves

CATALOG NO.	DESCRIPTION
MLS-PMT003	Triple Magnifier-Plastic. Handy molded plastic magnifier has 2X, 6X, and 8X magnification lenses. Overall length is 100mm (4").



CATALOG NO.	DESCRIPTION
MLS-PMD002	Clear Plastic Magnifier. Clear plastic magnifier is 110 mm (4.50") in length and features dual lenses of 3x and 6x magnification. 3x lens is 50 mm in diameter and 6x lens is 20mm in diameter.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-MPS038	Single Folding Magnifiers, Large, 5x
MLS-MPD038	Double Folding Magnifiers, Large, 10x
MLS-MPT038	Triple Folding Magnifiers, Large, 15x

CATALOG NO.	DESCRIPTION
MLS-MFG050	Hand Lens-Metal Frame, 50mm
MLS-MFG065	Hand Lens-Metal Frame, 65mm
MLS-MFG075	Hand Lens-Metal Frame, 75mm
MLS-MFG100	Hand Lens-Metal Frame, 100mm



CATALOG NO.	DESCRIPTION
MLS-MFP050	Hand Lens-Plastic Frame, 50mm
MLS-MFP065	Hand Lens-Plastic Frame, 65mm
MLS-MFP075	Hand Lens-Plastic Frame, 75mm
MLS-MFP090	Hand Lens-Plastic Frame, 90mm
MLS-MFP100	Hand Lens-Plastic Frame, 100mm

CATALOG NO.	DESCRIPTION
MLS-MFC010	Pocket Loupe, 10x magnification.



CATALOG NO.	DESCRIPTION
MLS-SGHC04	Density Blocks, Hooked Set of 4

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	MLS-SET462	MLS-SET583	MLS-SET694	MLS-SET615
DESCRIPTION	Support Stand and Ring Sets made from cast iron.			
BASE SIZE	4" x 6"	5" x 8"	6" x 9"	6" x 11"
ROD LENGTH	18"	20"	24"	36"
ROD DIAMETER	5/16"	3/8"	1/2"	1/2"
RING SIZE	3", 4"	3", 4", 5"	3", 4", 5", 6"	2", 3", 4", 5", 6"

CATALOG NO.	MLS-SSB4X6-CI	MLS-SSB5X8-CI	MLS-SSB6X9-CI	MLS-SSB6X11-CI
DESCRIPTION	Support Stands with Rods, Cast Iron			
BASE SIZE	4" x 6"	5" x 8"	6" x 9"	6" x 11"
ROD LENGTH	18"	20"	24"	36"
ROD DIAMETER	5/16"	3/8"	1/2"	1/2"



CATALOG NO.	DESCRIPTION
MLS-SGVC04	Equal Volume Metal Cylinders, Set of 4. Set of four includes aluminum, copper, steel, and brass cylinders. Includes a storage box and activity guide.

CATALOG NO.	DESCRIPTION
MLS-SGVS04	Equal Volume Metal Squares, Set of 4. Set of four includes aluminum, copper, steel, and brass cylinders. Includes a storage box and activity guide.



CATALOG NO.	DESCRIPTION
MLS-SGMC05	Equal Mass Metal Cylinders, Set of 5. Materials included are aluminum, zinc, tin, steel and copper. Includes a storage box and activity guide.

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-WSST09	Hooked Weight Sets, Stainless Steel. Set of 9 includes one each of 5g, 10g, 50g, 100g, and 500g, plus two each of 20g and 200g weights.
MLS-WSST10	Hooked Weight Sets, Stainless Steel. Set of 10 includes one each of 5g, 10g, 50g, 100g, 500g and 1000g, plus two each of 20g and 200g weights.

CATALOG NO.	DESCRIPTION
MLS-WHB005	Individual Hooked Brass Mass; 5G
MLS-WHB010	Individual Hooked Brass Mass; 10G
MLS-WHB020	Individual Hooked Brass Mass; 10G
MLS-WHB050	Individual Hooked Brass Mass; 50G
MLS-WHB100	Individual Hooked Brass Mass; 100G
MLS-WHB200	Individual Hooked Brass Mass; 200G
MLS-WHB500	Individual Hooked Brass Mass; 500G
MLS-WHB1000	Individual Hooked Brass Mass; 1000G



CATALOG NO.	DESCRIPTION
MLS-WHSBE9	Hooked Weight Set of 9, Black Enamel. Set of nine weights consisting of one each 10g, 50g, 100g, 500g, 1000g, and two each 20g and 200g weights.

CATALOG NO.	DESCRIPTION
MLS-WHSET10	Basic Hooked Weight Set of 10, Stainless Steel. Set of ten hooked weights are made of high quality stainless steel. Includes one of 200g, two of 100g, three of 50g, and four of 20g weights.



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-WHSET9	Hooked, Double ended Brass Weight. Set of nine weights consisting of one each 10g, 50g, 100g, 500g, 1,000g, and two each 20g and 200g weights.

CATALOG NO.	DESCRIPTION
MLS-WHST10	Deluxe Slotted Weight Set of 10, Brass. Includes rack and ten brass weights, 1 x 500g, 5 x 100g, 1 x 50g, 2 x 20g, 1 x 10g.



CATALOG NO.	DESCRIPTION
MLS-WHST12	Deluxe Slotted Weight Set of 12, Brass. Includes rack and twelve brass weights, 500g, 2 x 200g, 100g, 50g, 2 x 20g, 10g, 5g, 2 x 2g, 1g.

CATALOG NO.	DESCRIPTION
MLS-WHST13	Slotted Weight Set of 13, Brass. Includes hanger, case, and thirteen brass weights, 1 x 50g, 9 x 20g, 1 x 10g, 2 x 5g. Total mass is 250 grams.



CATALOG NO.	DESCRIPTION
MLS-WGCC44	Wire Gauze with Ceramic Center; 4"x4", Pk/10
MLS-WGCC55	Wire Gauze with Ceramic Center; 5"x5", Pk/10
MLS-WGCC66	Wire Gauze with Ceramic Center; 6"x6", Pk/10

CATALOG NO.	DESCRIPTION
MLS-WEC026	Connecting Wire, Enameled Copper Magnet; 26 gauge; 277'
MLS-WEC028	Connecting Wire, Enameled Copper Magnet; 28 gauge; 407'



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-3408	Clamp Holder Universal (Square Bosshead)

CATALOG NO.	DESCRIPTION
MLS-VCB001	Vernier Calipers, Stainless Steel



CATALOG NO.	DESCRIPTION
MLS-VCDT06	Dial Vernier Calipers

CATALOG NO.	DESCRIPTION
MLS-VCP001	Vernier Calipers, Polypropylene



CATALOG NO.	DESCRIPTION
MLS-WHXP05	Hexagonal Mass Set of 5, Plated. Includes one each of 100g, 500g, and 1000g, and two each of 200g masses.

CATALOG NO.	DESCRIPTION
MLS-WHX200	Hexagonal Mass, 200g
MLS-WHX500	Hexagonal Mass, 500g
MLS-WHX1000	Hexagonal Mass, 1000g



PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
MLS-HRDKIT3	Chemistry Labware Kit

CONTENTS	QUANTITY
Dropping Bottle, 30ml, Barnes	2
Dropping Pipettes, Plastic	6
Beaker, Glass, 50ml	1
Beaker, Glass, 100ml	1
Beaker, Glass, 250ml	1
Evaporating Dish, Porcelain, 70mm	1
Stirring Rods, Glass, 8"	6
Erlenmeyer Flask, Glass, 125ml	1
Erlenmeyer Flask, Glass, 250ml	1
Test Tubes, Glass, 18 x 150mm	6
Cylinder, Plastic (PP), 10ml	1
Cylinder, Plastic (PP), 50ml	1
Cylinder, Plastic (PP), 100ml	1

CONTENTS	QUANTITY
Forceps, Stainless, 5" long	1
Crucible Tongs, Steel, 9" long	1
Funnel, Glass, Long Stem, 90mm	1
Funnel, Plastic, 75mm dia.	1
Test Tube Clamp, Brass	1
Wire Gauze, 5" x 5",	1
Crucible, Porcelain, 25ml with Lid	1
Watch Glass, 75mm dia	1
Plastic Scoop	1
Spatula, stainless, 3" long Blade	1
Rubber Stoppers, Solid	3
Rubber Stoppers, One-hole	3
Rubber Stoppers, Two-hole	3

PHYSICS AND CHEMISTRY LAB SUPPLIES



CATALOG NO.	DESCRIPTION
HRDKIT2-C	Chemistry Hardware Assortment

CONTENTS		
Beaker Tongs	Clamp Holder, Right Angle	Support Stand with Rod
Burette Clamp	Clay Pipe Triangle, 2"	Test Tube Rack, Plastic
Burner Tripod, 3" ring, 9" legs	Crucible Tongs, Nickel-plated Steel	Test Tube Clamp
Clamp, Mohr's Pinchcock	Ring Support with Clamp, 4"	Wire Gauze with Ceramic Center

CATALOG NO.	DESCRIPTION
MLS-DSET14	Dissecting Instruments, Deluxe Set of 14



CONTENTS	
Scale, 6"	Scissors, Blunt/Sharp, 5"
Forceps, Blunt, 5"	Scissors, Sharp, 4"
Forceps, Sharp, Straight, 4"	Scalpel Handle, #4 w/ Blade
Needle, Arrow	Teasing Needle w/ Plastic Handle
Needle, Sharp, Straight	Spatula
Mall & Probe Seeker	Needle, Angular, Blunt
Scalpel w/ fixed blade, 5"	Dropper