

Whirl-pak Sampling Products



WHIRL-PAK BAGS

Whirl-Pak® sample bags are strong, sterile, and leakproof. Made of polyethylene, they are designed to carry solids or liquids to a laboratory for testing. Whirl-Pak™ Standard Sample Bags feature puncture proof tabs for protection from damage due to wire-end protrusion and leak-proof closures. Ideal for sampling and product analysis in food and beverage, dairy, water, medical, veterinary, environmental, soil and industrial markets.

Featured Products

- **Thio Bags**
- **Stand up Bags**
- **Filter Bags**
- **Write on Bags**

Thio bags: Whirl-Pak® Thio Sterilized Sampling Bags provide a sustainable, secure, contaminant free flexible container that ensures dependable analysis results. The Thio bag is a proven, economical and efficient way to collect, store and transport samples while reducing costs, space and waste to meet your needs.

Stand up bags: Whirl-Pak® Undosed Stand Up Sterilized Sampling Bags provide a secure, contaminant free flexible container that ensures dependable analysis results and meets your sustainability initiatives. You can count on Whirl-Pak® products as a proven, economical and efficient way to collect, store and transport samples.

Filter bags: Whirl-Pak® Filter Sterilized Sampling Bags provide a secure, contaminant free flexible container that ensures clean pipetting for dependable results. You can count on Whirl-Pak® products as a proven, economical and efficient way to collect, store and transport samples.

Write-on bags: Whirl-Pak® Write-On Sterilized Sampling Bags provide a secure, contaminant free flexible container that ensures dependable analysis results. You can count on Whirl-Pak® products as a proven, economical and efficient way to collect, store and transport samples.

WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01195	Whirl-Pak Write-On Bags - 55 oz. (1627) Box of 500 Yellow Tape



CATALOG NO.	DESCRIPTION
MLS-B01196	Whirl-Pak Write-On Bags - 24 oz. (710 ml) Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01067	Whirl-Pak Write-On Bags - 1 oz. (29 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01065	Whirl-Pak Write-On Bags - 18 oz. (532 ml) - Box of 500 - Yellow Tape

CATALOG NO.	DESCRIPTION
MLS-B01064	Whirl-Pak Write-On Bags - 18 oz. (532 ml) - Box of 500 - Yellow Tape



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01063	Whirl-Pak Bags - 24 oz. (710 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01062	Whirl-Pak Write-On Bags - 4 oz. (118 ml) - Box of 500 - Yellow Tape

CATALOG NO.	DESCRIPTION
MLS-B01040	Whirl-Pak Thio-Bags - 3.38 oz. (100 ml) - Box of 100



CATALOG NO.	DESCRIPTION
MLS-B01027	Whirl-Pak Bags - 42 oz. (1242 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01020	Whirl-Pak Bags - 24 oz. - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01018	Whirl-Pak Bags - 13 oz. (384 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01009	Whirl-Pak Bags - 2 oz. (58 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B00994	Whirl-Pak Bags - 36 oz. (1065 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B00992	Whirl-Pak Bags - 7 oz. (207 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B00990	Whirl-Pak Bags - 27 oz. (798 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01239	Whirl-Pak Homogenizer Blender Round-Bottom Bags - 52 oz. (1538) Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01245	Whirl-Pak Speci-Sponge Environmental Surface Sampling Bags - 18 oz. (532 ml) - Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01254	Whirl-Pak Thio-Bags - 10.1 oz. (300 ml) - Box of 100



CATALOG NO.	DESCRIPTION
MLS-B01297	Whirl-Pak Write-On Bags - 24 oz. (710 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01299	Whirl-Pak Speci-Sponge Environmental Surface Sampling Bags - 55 oz. (1627 ml) - Box of 100



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01323	Whirl-Pak Bags – 69 oz. (2041) – Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01339	Whirl-Pak Flat Wire Bags with Write-On Strip – 4 oz. (118 ml) – Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01341	Whirl-Pak Flat Wire Bags with Write-On Strip – 18 oz. (532 ml) Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01348	Whirl-Pak Homogenizer Blender Filter Bags – 24 oz. (710 ml) Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01350	Whirl-Pak Scoop Bags – 18 oz. (532 ml) – Box of 50



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01350B	Whirl-Pak Bags - 69 oz. (204l) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01364	Whirl-Pak Stand-Up Bags - 4 oz. (118 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01365	Whirl-Pak Stand-Up Bags - 18 oz. (532 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01385	Whirl-Pak Homogenizer Blender Filter Bags - 7 oz. (207 ml) - Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01392	Whirl-Pak Speci-Sponge Environmental Surface Sampling Bags with Sterile Glove - 18 oz. (532 ml) - Box of 100



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01401	Whirl-Pak Stand-Up Bags - 24 oz. (710 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01402	Whirl-Pak Stand-Up Thio-Bags - 3.4 oz. (100 ml) - Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01403	Whirl-Pak Stand-Up Thio-Bags - 10.1 oz. (300 ml) - Box of 100



CATALOG NO.	DESCRIPTION
MLS-B01416	Whirl-Pak Homogenizer Blender Filter Bags - 69 oz. (2041 ml) - Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01418	Whirl-Pak Light Sensitive Black Bags - 18 oz. (532 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01421	Whirl-Pak Plain Blender Bags - 7-12 x 12 (19 x 30 cm) - Box of 1000



CATALOG NO.	DESCRIPTION
MLS-B01422	Whirl-Pak Hydrated Spec-Sponge Bags - 18 oz. (532 ml) - Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01423	Whirl-Pak Hydrated Spec-Sponge Bags with Sterile Glove - 18 oz. (532 ml) - Box of 100



CATALOG NO.	DESCRIPTION
MLS-B01445	Whirl-Pak Bags - 92 oz. (2721 ml) - Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01446	Whirl-Pak Bags - 123 oz. (3637 ml) - Box of 250



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01447	Whirl-Pak Bags - 184 oz. (5441 ml) - Box of 100



CATALOG NO.	DESCRIPTION
MLS-B01449	Whirl-Pak Stand-Up Bags - 36 oz. (1065 ml) - Box of 250



CATALOG NO.	DESCRIPTION
MLS-B01450	Whirl-Pak Stand-Up Bags - 42 oz. (1242 ml) - Box of 250



CATALOG NO.	DESCRIPTION
MLS-B01451	Whirl-Pak Stand-Up Bags - 69 oz. (2041 ml) - Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01472	Whirl-Pak Light SensitiveBlack Bags - 4 oz. (118 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01474	Whirl-Pak Plain Blender Bags - 15 x 20 (37.5 x 50 cm) - Box of 250



CATALOG NO.	DESCRIPTION
MLS-B01488	Whirl-Pak Homogenizer Blender Filter Bags - 92 oz. (2721 ml) Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01489	Whirl-Pak Write-On Bags - 7 oz. (207 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01490	Whirl-Pak Write-On Bags - 13 oz. (384 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01491	Whirl-Pak Write-On Bags - 27 oz. (798 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01515	Whirl-Pak Write-On Bags - 69 oz. (2041 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01525	Whirl-Pak Filter Bags - 138 oz. (4080 ml) - Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01532	Whirl-Pak Bags - 55 oz. (1627 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01542	Whirl-Pak Large Stand-Up Bags - 138 oz. (4080 ml) - Box of 100

CATALOG NO.	DESCRIPTION
MLS-B01545	Whirl-Pak Detectable Bags - 18 oz. (532 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01546	Whirl-Pak Detectable Bags - 42 oz. (1242 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01558	Whirl-Pak Light Sensitive Black Bags - 55 oz. (1627 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B01590	Whirl-Pak Hydrated PolySponge Bags Without Glove - 100 Sponges



CATALOG NO.	DESCRIPTION
MLS-B01591	Whirl-Pak Hydrated PolySponge Bags With Sterile Glove - 100 Sponges

CATALOG NO.	DESCRIPTION
MLS-B01601	Whirl-Pak Thio-Bag 500ml Stand-Up



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01603	Whirl-Pak Thio-Bag 1L Stand-Up



CATALOG NO.	DESCRIPTION
MLS-B01604	Whirl-Pak Thio-Bag 2L Stand-Up

CATALOG NO.	DESCRIPTION
MLS-B00679	Whirl-Pak Bags - 4 oz. (118 ml) - Box of 500



CATALOG NO.	DESCRIPTION
MLS-B00736	Whirl-Pak Bags - 18 oz. (532 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS- B01253	Whirl-Pak® Trans-Pak® Non-Sterile Specimen Bags - 24 oz. (710 ml) - Box of 500



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01547	Whirl-Pak Filter Bags Without Tape and Wire - 55 oz. (1627 ml) - Box of 250



CATALOG NO.	DESCRIPTION
MLS-B01318	Whirl-Pak Homogenizer Blender Filter Bags - 55 oz. (1627 ml) - Box of 250

CATALOG NO.	DESCRIPTION
MLS-B01592	Whirl-Pak® Hydrated PolyProbe™ - 100 Probes



CATALOG NO.	DESCRIPTION
MLS-B01255	Whirl-Pak® Pocket Bags - 24 oz. (710 ml) - Box of 500

CATALOG NO.	DESCRIPTION
MLS- B01318	Whirl-Pak® Homogenizer Blender Filter Bags - 55 oz. (1,627 ml) - Box of 250



WHIRL-PAK BAGS

CATALOG NO.	DESCRIPTION
MLS-B01319	Nasco Poultry Rinse Bags - 15 in. W x 20 in. L (37.5 x 50 cm) - Box of 250



CATALOG NO.	DESCRIPTION
MLS-B01324	Nasco Meat/Poultry Carcass Sampling Kit



CATALOG NO.	DESCRIPTION
MLS-B01475	Whirl-Pak® Sterile Dry Sponge Probe



CATALOG NO.	DESCRIPTION
MLS-B01478	Whirl-Pak® Spoon Bag - 18 oz. (532 ml) - Box of 100

WHIRL-PAK RACKS




CATALOG NO.	DESCRIPTION
MLS-B00677	(DNR), Rack, 21 Compartment W/Double Grid

CATALOG NO.	DESCRIPTION
MLS-B00678	(DNR), Rack, 15 Compartment W/Double Grid




CATALOG NO.	DESCRIPTION
MLS-B00750	(DNR), Rack, 30 Compartment W/Double Grid

CATALOG NO.	DESCRIPTION
MLS-B00751	Rack 12, Compartment W/Double Grid



CATALOG NO.	DESCRIPTION
MLS-B00772	(DNR), Rack, 18 Compartment W/Double Grid

CATALOG NO.	DESCRIPTION
MLS-B01014	Rack, 20Compartment



WHIRL-PAK RACKS



CATALOG NO.	DESCRIPTION
MLS-B01023	(DNR), Rack, 18 Compartment Single Grid

CATALOG NO.	DESCRIPTION
MLS-B01048	(DNR), Rack, 15 Compartment W/Double Grid



CATALOG NO.	DESCRIPTION
MLS-B01109	(DNR), Rack, 15 Compartment Double Grid

CATALOG NO.	DESCRIPTION
MLS-B01533	Rack for Bag B01445 with Handle



CATALOG NO.	DESCRIPTION
MLS-B01534	(DNR) Rack for Bag B01446 with Handle

CATALOG NO.	DESCRIPTION
MLS-B01535	(DNR) Rack for Bag B01447 with Handle



SLUDGE JUDGE



CATALOG NO.	DESCRIPTION
MLS-C09247	Sludge Judge® - 15 ft. (4.60 m)



CATALOG NO.	DESCRIPTION
MLS-C09250	Sludge Judge® Bottom Section with Valve 5 ft. (1.53 m)



CATALOG NO.	DESCRIPTION
MLS-C09251	Sludge Judge® Extension Section - 5 ft. (1.53 m)



CATALOG NO.	DESCRIPTION
MLS-C09252	Sludge Judge® Top Section with Rope 5 ft. (1.53 m)



CATALOG NO.	DESCRIPTION
MLS-B01355	Sludge Judge® II - Large Capacity

CATALOG NO.	DESCRIPTION
MLS-B01356	Sludge Judge® II Bottom Section with Valve - 5 ft. (1.53 m)



SLUDGE JUDGE



CATALOG NO.	DESCRIPTION
MLS- B01357	Sludge Judge® II Extension Section - 5 ft. (1.53 m)

CATALOG NO.	DESCRIPTION
MLS-B01358	Sludge Judge® II Top Section with Rope - 5 ft. (1.53 m)



CATALOG NO.	DESCRIPTION
MLS-B01386	Sludge Judge® Ultra - Super Strong

CATALOG NO.	DESCRIPTION
MLS-B01387	Sludge Judge® Ultra Top Section with Rope - 5 ft. (1.53 m)



CATALOG NO.	DESCRIPTION
MLS-B01388	Sludge Judge® Ultra Top Section with Rope - 5 ft. (1.53 m)

CATALOG NO.	DESCRIPTION
MLS-B01389	Sludge Judge® Ultra Top Section with Rope - 5 ft. (1.53 m)



SAMPLER CARE TOOLS



CATALOG NO.	DESCRIPTION
MLS- B01291	Cleaning Brush - 6 ft. (1.83 m) L. 1 in. (2.54 cm) Diameter Bristle
MLS- B01361	Cleaning Brush 6Ft For B1355

CATALOG NO.	DESCRIPTION
MLS- C09249	Case for Sampler



COLLECTION POLES



CATALOG NO.	DESCRIPTION
MLS- B01229	Nasco Sludge Nabber

CATALOG NO.	DESCRIPTION
MLS-B01366	Nasco Swing Sampler - 24 ft. (7.32 m)



CATALOG NO.	DESCRIPTION
MLS-B01310	Nasco Swing Sampler - 12 ft. (3.66 m)

CATALOG NO.	DESCRIPTION
MLS-B01354	Replacement Bottle Holder for Nasco Swing Samplers



ACCESSORIES



CATALOG NO.	DESCRIPTION
MLS- B01595	Locking Pipe Closures - Pack of 50 - Yellow and White
MLS-B01550	Locking Pipe Closures - Pack of 50



CATALOG NO.	DESCRIPTION
MLS-B01523	Whirl-Pak® Pipette Filter Sleeve - Box of 500

CATALOG NO.	DESCRIPTION
MLS-B01559	Replacement Sponge Holder for Nasco Sponge Pole



CATALOG NO.	DESCRIPTION
MLS-B01593	Whirl-Pak Probe Adapter